

TMEM9 antibody



Catalog Number: 116199

Product name

TMEM9 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TMEM9 Rabbit Polyclonal antibody. Positive WB detected in mouse brain tissue, Transfected HEK-293 cells. Observed molecular weight by Western-blot: 30 kDa

Preparation

This antibody was obtained by immunization of TMEM9 recombinant protein (Accession Number: NM_016456). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB

Dilutions


Recommended Dilution:

WB: 1:200-1:1000

Validations



Transfected HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:116199(TMEM9 Antibody) at dilution of 1:700



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:116199(TMEM9 Antibody) at dilution of 1:300