

# Anti-PRMT6 (2C3) Mouse antibody



Catalog Number: 168199

## Product name

Anti-PRMT6 (2C3) Mouse antibody

## Specificity

Human

## Antibody description

PRMT6 (2C3) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human PRMT6 protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal in buffer containing 0.1M Tris-Glycine (pH 7.4 150 mM NaCl) with 0.02% sodium azide 50% glycerol

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG1

## Applications


WB, IHC-P, ICC/IF

## Dilutions


WB: 1/1000

;ICC: 1/300

## Validations




Immunohistochemistry stain of paraffin-embedded human breast cancer using PRMT6 mouse mAb



(1:200).

Western blot detection of PRMT6 in HeLa and 293T cell lysates using PRMT6 mouse mAb (1:1000 diluted). Predicted band size: 42KDa. Observed band size: 42KDa.



Immunocytochemistry stain of HeLa using PRMT6 mouse mAb (1:300).