

# Human PARVA/MXRA2 (His tag) recombinant protein



Catalog Number: 518199

## General Information

### Gene Name Synonym

Alpha-Parvin, Actopaxin, CH-ILKBP, Calponin-Like Integrin-Linked Kinase-Binding Protein, Matrix-Remodeling-Associated Protein 2, PARVA, MXRA2

### Protein Construction

A DNA sequence encoding the Human PARVA/MXRA2 Met1-Glu372 is expressed with a 6His tag at the C-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 50mM Tris, 150mM NaCl, 40% Glycerol, pH 7.4.

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Human PARVA/MXRA2 (His tag) recombinant protein