

Human NPM1/Nucleophosmin (His Tag) recombinant protein



Catalog Number: 500599

General Information

Gene Name Synonym

Nucleolar phosphoprotein B23; Nucleolar protein NO38; Numatrin

Protein Construction

A DNA sequence encoding the human NPM1 isoform 1 (P06748-1) (Glu2-Leu294) was expressed with a polyhistidine tag at the N-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

> 85 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant human NPM1 consists of 304 amino acids and has a calculated molecular mass of 34 kDa.

Formulation

Lyophilized from sterile 30 mM Hepes, 2 mM EDTA, 0.001 % Tween, 15 % glycerol, pH 7.0. 1. 5 % trehalose and mannitol are added as protectants before lyophilization. 2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human NPM1 / Nucleophosmin (aa 2-294) Protein (His Tag)