

# Human PILR-alpha/PILRA recombinant protein



Catalog Number: 502599

## General Information

### Gene Name Synonym

Cell surface receptor FDF03; Inhibitory receptor PILR-alpha

### Protein Construction

A DNA sequence encoding the human PILRA (AAH17812.1)(Met1-Asp226) was expressed with five amino acids (DDDDK) at the C-terminus.

### Organism

Human

### Expression Host

Human Cells

## QC Testing

### Purity

> 85 % as determined by SDS-PAGE

### Endotoxin

< 1.0 EU per  $\mu\text{g}$  of the protein as determined by the LAL method

### Stability

Samples are stable for up to twelve months from date of receipt at  $-70^{\circ}\text{C}$

### Predicted N terminal

Gln 20

### Molecular Mass

The recombinant human PILRA comprises 213 amino acids and has a predicted molecular mass of 24.1 kDa. The apparent molecular mass of the protein is approximately 33 kDa in SDS-PAGE under reducing conditions.

### Formulation

Lyophilized from sterile PBS, pH 7.4.

1. 5 % trehalose and mannitol are added as protectants before lyophilization.
2. Please contact us for any concerns or special requirements.

## Usage Guide


### Storage

Store it under sterile conditions at  $-20^{\circ}\text{C}$  to  $-80^{\circ}\text{C}$ . It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

### Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

### SDS-PAGE



Human PILR-alpha / PILRA Protein SDS-PAGE