

SNRPD2 antibody



Catalog Number: 107599

Product name

SNRPD2 antibody

Specificity

Human; other species not tested.

Antibody description

SNRPD2 Mouse Monoclonal antibody. Positive WB detected in Hela cells, MCF7 cells. Observed molecular weight by Western-blot: 16kDa

Preparation

This antibody was obtained by immunization of SNRPD2 recombinant protein (Accession Number: NM_001384647). Purification method: Protein A purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Monoclonal

Ig Type

Mouse IgG2b

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:107599(SNRPD2 antibody) at dilution of 1:1000