

# BRUNOL5 antibody



Catalog Number: 108599

## Product name

BRUNOL5 antibody

## Specificity

human, mouse, rat; other species not tested.

## Antibody description

BRUNOL5 Rabbit Polyclonal antibody. Positive IHC detected in human gliomas. Positive WB detected in mouse skeletal muscle tissue, mouse testis tissue. Observed molecular weight by Western-blot: 52kd

## Preparation

This antibody was obtained by immunization of BRUNOL5 recombinant protein (Accession Number: XM\_006722832). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


## Validations



mouse skeletal muscle tissue were subjected to SDS PAGE followed by western blot with Catalog No:117246(BRUNOL5 antibody) at dilution of



1:1500  
Immunohistochemical of paraffin-embedded human gliomas using Catalog No:117246(BRUNOL5 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human gliomas using Catalog No:117246(BRUNOL5 antibody) at dilution of 1:50

# BRUNOL5 antibody



Catalog Number: 108599

---

(under 40x lens)