

PSEN1 antibody

Catalog Number: 114299

Product name

PSEN1 antibody

Specificity

Human; other species not tested.

Antibody description

PSEN1 Rabbit Polyclonal antibody. Positive IF detected in HEK-293 cells.

Preparation

This antibody was obtained by immunization of PSEN1 recombinant protein (Accession Number: NM_007318). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, IF

Dilutions

Recommended Dilution:

IF: 1:20-1:200

Validations

Immunofluorescent analysis of HEK-293 cells using Catalog No:114299(PSEN1 Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)