

SLC12A4 antibody



Catalog Number: 115299

Product name

SLC12A4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SLC12A4 Rabbit Polyclonal antibody. Positive IP detected in HeLa cells. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 121 kDa

Preparation

This antibody was obtained by immunization of SLC12A4 recombinant protein (Accession Number: NM_005072). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IP: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115299(SLC12A4 antibody) at dilution of 1:400



IP Result of anti-SLC12A4 (IP:Catalog No:115299, 4ug; Detection:Catalog No:115299 1:500) with HeLa cells lysate 3200ug.