

# Mouse anti V5-Tag mAb



Catalog Number: 166299

## Product name

Mouse anti V5-Tag mAb

## Antibody description

Monoclonal antibody to Mouse anti V5-Tag

## Preparation

Antigen: A synthetic peptide of human V5-Tag

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Monoclonal

## Ig Type

Mouse IgG

## Applications

WB, IF, IP


## Dilutions

WB 1:1000 - 1:5000

IF 1:50 - 1:200

IP 1:50 - 1:200

## Validations



Western blot - Mouse anti V5-Tag mAb

Western blot analysis of recombinant V5 protein with Mouse anti V5-Tag mAb at dilution of 1:2000. Secondary antibody: HRP Goat Anti-Mouse IgG (H+L) at 1:10000 dilution. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 1s.