

General Information

Gene Name Synonym

FIL1 eta; Interleukin-1 eta; Interleukin-1 family member 8; Interleukin-1 homolog 2

Protein Construction

A DNA sequence encoding the mature form of human IL36B isoform 2 (Q9NZH7-2) (Met 1-Glu 157) was expressed.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

> 97 % as determined by SDS-PAGE

Endotoxin

Please contact us for more information.

Stability

Samples are stable for up to twelve months from date of receipt at -70°C

Predicted N terminal

Met

Molecular Mass

The recombinant human IL36B consists of 157 amino acids and has a predicted molecular mass of 17.7 kDa. It migrates as an approximately 18 kDa band in SDS-PAGE under reducing conditions.

Formulation

Lyophilized from sterile PBS, pH 7.4

1. 5 % trehalose and mannitol are added as protectants before lyophilization.

2. Please contact us for any concerns or special requirements.

Usage Guide


Storage

Store it under sterile conditions at -20°C to -80°C. It is recommended that the protein be aliquoted for optimal storage. Avoid repeated freeze-thaw cycles.

Reconstitution

Adding sterile water, prepare a stock solution of 0.25 mg/ml. Concentration is measured by UV-Vis.

SDS-PAGE



Human IL36B / IL1F8 Protein SDS-PAGE