

UVRAG antibody



Catalog Number: 116699

Product name

UVRAG antibody

IP: 1:200-1:2000

IF: 1:10-1:100

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

UVRAG Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive IP detected in mouse liver tissue. Positive WB detected in HepG2 cells, HEK-293 cells, human brain tissue, mouse liver tissue. Observed molecular weight by Western-blot: 75-85 kDa

Preparation

This antibody was obtained by immunization of UVRAG recombinant protein (Accession Number: NM_003369). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB, IF, IP

Dilutions


Recommended Dilution:

WB: 1:500-1:5000


Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:116699(UVRAG antibody) at dilution of 1:500



Immunofluorescent analysis of HepG2 cells, using UVRAG antibody Catalog No:116699 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit



IP Result of anti-UVRAG (IP:Catalog No:116699,

UVRAG antibody



Catalog Number: 116699

4ug; Detection:Catalog No:116699 1:600) with mouse liver tissue lysate 4000ug.