

ACSM2A antibody



Catalog Number: 107699

Product name

Polyclonal

ACSM2A antibody

Ig Type

Specificity

Rabbit IgG

Human; other species not tested.

Applications

Antibody description

ELISA, WB

ACSM2A Rabbit Polyclonal antibody. Positive WB detected in human brain tissue. Observed molecular weight by Western-blot: 64 kDa

Dilutions


Recommended Dilution:

WB: 1:500-1:5000

Preparation

This antibody was obtained by immunization of ACSM2A recombinant protein (Accession Number: NM_001010845). Purification method: Antigen affinity purified.

Validations



Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

human brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:107699(ACSM2A antibody) at dilution of 1:1000