

Anti-p53 Rabbit antibody



Catalog Number: 168699

Product name

Anti-p53 Rabbit antibody

Specificity

Human

Antibody description

p53 Rabbit monoclonal antibody

Preparation

Antigen: Recombinant protein of human p53

Formulation

50nM Tris-Glycine(pH 7.4), 0.15M NaCl,
40%Glycerol, 0.01% sodium azide and 0.05% BSA.

Storage

Store at 4°C short term. Store at -20°C long term.
Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications


WB, ICC/IF

Dilutions

WB: 1/1000

ICC/IF: 1/50

Validations



Western blot detection of p53 in Jurkat cell lysates using p53 Rabbit mAb(1:500 diluted).Predicted band size:44kDa.Observed band size:53kDa.