

Human CYGB/STAP (His tag) recombinant protein



Catalog Number: 518699

General Information

Gene Name Synonym

Cytoglobin, Histoglobin, HGb, Stellate Cell Activation-Associated Protein, CYGB, STAP

Protein Construction

A DNA sequence encoding the Human CYGB/STAP Met1-Pro190 is expressed with a 6His tag at the C-terminus.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL

method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM TrisHCl, 10% Glycerol, pH 8.0 .

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human CYGB/STAP (His tag) recombinant protein