

Anti-OTUB2 antibody



Catalog Number: 100399

Product name

Anti-OTUB2 antibody

Immunogen

[Human OTUB2 \(His Tag\) recombinant protein](#)

Specificity

Human OTUB2

Antibody description

Rabbit monoclonal to OTUB2

Preparation

This antibody was obtained from a rabbit immunized with purified, recombinant Human OTUB2 (rh OTUB2; Q96DC9-1; Met1-His234).

Formulation

0.2 μ m filtered solution in PBS

Storage

This antibody can be stored at 2°C-8°C for one month without detectable loss of activity. Antibody products are stable for twelve months from date of receipt when stored at -20°C to -80°C.

Preservative-Free.

Sodium azide is recommended to avoid contamination (final concentration 0.05%-0.1%). It is toxic to cells and should be disposed of properly. Avoid repeated freeze-thaw cycles.

Clonality

Monoclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF, ICC/IF


Dilutions

ELISA: 0.1-0.2 μ g/ml

This antibody can be used at 0.1-0.2 μ g/ml with the appropriate secondary reagents to detect Human OTUB2. The detection limit for Human OTUB2 is approximately 0.00245 ng/well.

ICC/IF: 10-25 μ g/mL

Validations



OTUB2 Antibody, Rabbit MAb

Immunofluorescence staining of Human OTUB2 in A549 cells. Cells were fixed with 4% PFA, permeabilized with 0.3% Triton X-100 in PBS, blocked with 10% serum, and incubated with rabbit anti-Human OTUB2 monoclonal antibody (15 μ g/ml) at 4°C overnight. Then cells were stained with the Alexa Fluor® 549-conjugated Goat Anti-rabbit IgG secondary antibody (red). Positive staining was localized to cytoplasm and nucleus.