

Anti-B4GALT5 antibody



Catalog Number: 176399

Product name

Anti-B4GALT5 antibody

Ig Type

Rabbit IgG

Specificity

Human, Mouse

Applications

WB

Antibody description

Rabbit polyclonal antibody to B4GALT5

Dilutions

WB: 1:500-1:2000

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human b4galt5. aa range:321-370

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Validations




Fig1:; Western Blot analysis of various cells using β -1,4-Gal-T5 Polyclonal Antibody diluted at 1:1000. Secondary antibody (catalog#: HA1001) was diluted at 1:20000