

Anti-B4GALT5 antibody

Catalog Number: 176399

Product name

Anti-B4GALT5 antibody

Specificity

Human, Mouse

Antibody description

Rabbit polyclonal antibody to B4GALT5

Preparation

This antigen of this antibody was the antiserum was produced against synthesized peptide derived from human b4galt5. aa range:321-370

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1:; Western Blot analysis of various cells using β -1,4-Gal-T5 Polyclonal Antibody diluted at 1:1000. Secondary antibody (catalog#: HA1001) was diluted at 1:20000