

Anti-IL-8 (6E7) Mouse antibody



Catalog Number: 168399

Product name

Anti-IL-8 (6E7) Mouse antibody

Specificity

Human

Antibody description

IL-8 (6E7) Mouse monoclonal antibody

Preparation

Antigen: Purified recombinant human IL8 protein fragments expressed in E.coli.

Formulation

Purified mouse monoclonal in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

Clonality

Monoclonal

Ig Type

Mouse IgG2b


Applications

WB, ELISA


Dilutions

WB: 1/500

Validations



Indirect ELISA assay for anti-IL8 mouse mAb. Antigen coating concentration: 8ug/ml.



Western blot detection of IL8 in CHO-K1 transfected by IL8-PDGFR fusion protein cell lysate using IL8 mouse mAb (1:500(A)-1:1000(B) diluted).