

Human ANXA5/ANX5/ENX2/PP4 recombinant protein



Catalog Number: 519399

General Information

Gene Name Synonym

Annexin A5, Anchorin CII, Annexin V, Annexin-5, Calphobindin I, CBP-I, Endonexin II, Lipocortin V, Placental Anticoagulant Protein 4, PP4, Placental Anticoagulant Protein I, PAP-I, Thromboplastin Inhibitor, Vascular Anticoagulant-Alpha, VAC-Alpha, ANXA5, ANX5, ENX2, PP4

Protein Construction

A DNA sequence encoding the Human ANXA5/ANX5/ENX2/PP4 Ala2-Asp320 is expressed.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM TrisHCl, 1mM DTT, 10% Glycerol, pH 8.0.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human ANXA5/ANX5/ENX2/PP4 recombinant protein