

RPAP3 antibody



Catalog Number: 114799

Product name

RPAP3 antibody

Specificity

Human; other species not tested.

Antibody description

RPAP3 Rabbit Polyclonal antibody. Positive WB detected in HEK-293 cells, HeLa cells, Jurkat cells. Observed molecular weight by Western-blot: 76 kDa

Preparation

This antibody was obtained by immunization of RPAP3 recombinant protein (Accession Number: NM_001146075). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, WB

Dilutions

Recommended Dilution:

WB: 1:500-1:5000

Validations



HEK-293 cells were subjected to SDS PAGE followed by western blot with Catalog No:114799(RPAP3 antibody) at dilution of 1:1000