

C9orf41 antibody



Catalog Number: 108799

Product name

C9orf41 antibody

Immunogen

[Human C9orf41 Recombinant protein \(His tag\)](#)

Specificity

human; other species not tested.

Antibody description

C9orf41 Rabbit Polyclonal antibody. Positive WB detected in HepG2 cells, A549 cells, HeLa cells, MCF7 cells.

Preparation

This antibody was obtained by immunization of C9orf41 recombinant protein (Accession Number: NM_152420). Purification method: Antigen Affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG


Applications

ELISA, WB

Dilutions

Recommended Dilution:
WB: 1:500-1:5000

Validations



HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:108799(C9orf41 antibody) at dilution of 1:1000