

Human HMOX1/HO/HO1 recombinant protein



Catalog Number: 518799

General Information

Gene Name Synonym

Heme Oxygenase 1, HO-1, HMOX1, HO, HO1

Protein Construction

A DNA sequence encoding the Human HMOX1/HO/HO1 Met1-Thr261 is expressed.

Organism

Human

Expression Host

E. coli

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM PB, 150mM NaCl, 1mM EDTA, pH 7.4.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Human HMOX1/HO/HO1 recombinant protein