

GNG4 antibody

Catalog Number: 111002

Product name

GNG4 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

GNG4 Rabbit Polyclonal antibody. Positive IHC detected in human pancreas cancer tissue. Positive WB detected in mouse brain tissue. Observed molecular weight by Western-blot: 8 kDa

Preparation

This antibody was obtained by immunization of GNG4 recombinant protein (Accession Number: NM_001098722). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC

Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IHC: 1:20-1:200


Validations



mouse brain tissue were subjected to SDS PAGE followed by western blot with Catalog No:111002(GNG4 antibody) at dilution of 1:400



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:111002(GNG4 antibody) at dilution of 1:100 (under 10x lens)



Immunohistochemical of paraffin-embedded human pancreas cancer using Catalog No:111002(GNG4 antibody) at dilution of 1:100

GNG4 antibody



Catalog Number: 111002

(under 40x lens)