

# TONSL Polyclonal Antibody



Catalog Number: 161402

## Product name

TONSL Polyclonal Antibody

## Specificity

Human, Mouse, Rat

## Antibody description

Polyclonal antibody to TONSL

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 160-390 of human TONSL (NP\_038460.4).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


WB, IHC

## Dilutions

WB 1:500 - 1:2000


IHC 1:100 - 1:200

## Validations



### Western blot - TONSL Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using TONSL antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 30s.



### Western blot - TONSL Polyclonal Antibody

Western blot analysis of extracts of mouse brain, using TONSL antibody at 1:3000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.