

# ACADSB Polyclonal Antibody



Catalog Number: 162402

---

## Product name

ACADSB Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to ACADSB

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 133-432 of human ACADSB (NP\_001600.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IF

## Dilutions

IF 1:50 - 1:200

## Validations



Immunofluorescence - ACADSB Polyclonal Antibody

Immunofluorescence analysis of HeLa cells using ACADSB antibody at dilution of 1:100 (40x lens). Blue: DAPI for nuclear staining.