

# LY6G6F antibody



Catalog Number: 112402

## Product name

LY6G6F antibody

## Specificity

Human; other species not tested.

## Antibody description

LY6G6F Rabbit Polyclonal antibody. Positive IF detected in HepG2 cells. Positive WB detected in HepG2 cells, HL-60 cells. Observed molecular weight by Western-blot: 60 kDa

## Preparation

This antibody was obtained by immunization of LY6G6F recombinant protein (Accession Number: NM\_001003693). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IF


## Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:20-1:200

## Validations




HepG2 cells were subjected to SDS PAGE followed by western blot with Catalog No:112402(LY6G6F



Antibody) at dilution of 1:1000

HL-60 cells were subjected to SDS PAGE followed by western blot with Catalog No:112402(LY6G6F Antibody) at dilution of 1:1000



Immunofluorescent analysis of HepG2 cells using Catalog No:112402(LY6G6F Antibody) at dilution of 1:50 and Alexa Fluor 488-conjugated AffiniPure Goat Anti-Rabbit IgG(H+L)