

# CADPS Polyclonal Antibody



Catalog Number: 163402

## Product name

CADPS Polyclonal Antibody

## Specificity

Human

## Antibody description

Polyclonal antibody to CADPS

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1020-1100 of human CADPS (NP\_899630.1).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

IHC

## Dilutions


IHC 1:50 - 1:200

## Validations



Immunohistochemistry - CADPS Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human breast cancer using CADPS antibody at dilution of 1:100 (40x lens).



Immunohistochemistry - CADPS Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using CADPS antibody at dilution of 1:100 (40x lens).