

UBE2W Polyclonal Antibody



Catalog Number: 164402

Product name

UBE2W Polyclonal Antibody

Specificity

Human, Mouse

Antibody description

Polyclonal antibody to UBE2W

Preparation

Antigen: A synthetic Peptide of human UBE2W

Formulation

PBS with 0.02% sodium azide, pH7.3.

Storage

Store at 4°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


WB, IHC

Dilutions

WB 1:500 - 1:1000

IHC 1:20 - 1:50

Validations



Western blot - UBE2W Polyclonal Antibody

Western blot analysis of extracts of MCF-7 cells, using UBE2W antibody .Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution.Lysates/proteins: 25ug per lane.Blocking buffer: 3% nonfat dry milk in TBST.