

# Anti-5 Lipoxygenase Rabbit antibody



Catalog Number: 167402

## Product name

Rabbit IgG

Anti-5 Lipoxygenase Rabbit antibody

## Applications

## Specificity

WB, IHC-P, ICC/IF, FC

Human, Mouse, Rat

## Dilutions

## Antibody description

WB: 1/500-1/2000

5 Lipoxygenase Rabbit polyclonal antibody

IHC: 1/50-1/200

## Preparation

ICC/IF: 1/50-1/200


Antigen: A synthesized peptide derived from human 5 Lipoxygenase

FC: 1/50

## Validations

## Formulation

Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.



## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Polyclonal

Western blot detection of 5 Lipoxygenase in Raji, MCF7 cell lysates using 5 Lipoxygenase Rabbit pAb (1:1000 diluted). Predicted band size: 78kDa. Observed band size: 78kDa.

## Ig Type