

# Anti-PRMT1 (4E9) Mouse antibody



Catalog Number: 168402

## Product name

Anti-PRMT1 (4E9) Mouse antibody

## Specificity

Human

## Antibody description

PRMT1 (4E9) Mouse monoclonal antibody

## Preparation

Antigen: Purified recombinant human PRMT1 protein fragments expressed in E.coli.

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.02% sodium azide and 50% glycerol.

## Storage

Store at 4°C short term. Store at -20°C long term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal

## Ig Type

Mouse IgG2b


## Applications

WB

## Dilutions

WB: 1/1000

## Validations



Western blot detection of PRMT1 in HeLa,A549 and SW480 cell lysates using PRMT1 mouse mAb (1:1000 diluted).Predicted band size:42KDa.Observed band size:42KDa.