

# MCM3AP Polyclonal Antibody



Catalog Number: 162802

## Product name

MCM3AP Polyclonal Antibody

## Specificity

Mouse

## Antibody description

Polyclonal antibody to MCM3AP

## Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 250-550 of human MCM3AP (NP\_003897.2).

## Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

## Storage

Store at -20°C. Avoid freeze / thaw cycles.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG


## Applications

WB

## Dilutions

WB 1:500 - 1:2000

## Validations



## Western blot - MCM3AP Polyclonal Antibody

Western blot analysis of extracts of mouse spleen, using MCM3AP antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Enhanced Kit. Exposure time: 60s.