

# Anti-Bak Rabbit antibody



Catalog Number: 166802

## Product name

Anti-Bak Rabbit antibody

IHC: 1/20

ICC/IF: 1/50

## Specificity

Human

FC: 1/20

IP: 1/20

## Antibody description

Bak Rabbit polyclonal antibody

## Validations

## Preparation

Antigen: A synthetic peptide of human Bak

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type


Rabbit IgG

## Applications


WB, IHC-P, ICC/IF, FC, IP

## Dilutions

WB: 1/1000



Western blot detection of Bak in K562, HeLa cell lysates using Bak Rabbit pAb(1:1000 diluted). Predicted band size: 23kDa. Observed band size: 23kDa.



Immunohistochemistry of Bak in paraffin-embedded Human colon cancer tissue using Bak Rabbit pAb at dilution 1/20