

Mouse Nt5e/Nt5/Nte (His tag) recombinant protein



Catalog Number: 518802

General Information

Gene Name Synonym

5'-nucleotidase, Ecto-5'-nucleotidase, CD73, 5'-NT

Protein Construction

A DNA sequence encoding the Mouse Nt5e/Nt5/Nte Trp29-Phe550 is expressed with a 6His tag at the C-terminus.

Organism

Mouse

Expression Host

Human Cells

QC Testing

Purity

>95% as determined by SDS-PAGE.

Endotoxin

<1.0 EU per μg protein as determined by the LAL method.

Stability

Store for up to 6 months at -20°C

Formulation


Supplied as a $0.2\ \mu\text{m}$ filtered solution of 20mM TrisHCl, 120mM NaCl, 4mM CaCl₂, 20% Glycerol, pH 7.5.

Usage Guide

Storage

Store at $<-20^{\circ}\text{C}$, minimize freeze-thaw cycles.

SDS-PAGE



Mouse Nt5e/Nt5/Nte (His tag) recombinant protein