

PTHLH Polyclonal Antibody



Catalog Number: 161102

Product name

PTHLH Polyclonal Antibody

WB 1:500 - 1:2000

IHC 1:50 - 1:200

Specificity

Human

Antibody description

Polyclonal antibody to PTHLH

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 74-130 of human PTHLH (NP_945316.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type


Rabbit IgG

Applications

WB, IHC


Dilutions

Validations



Immunohistochemistry - PTHLH Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human prostate using PTHLH antibody at dilution



of 1:100 (40x lens).

Immunohistochemistry - PTHLH Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human stomach using PTHLH antibody at dilution of 1:100 (40x lens).