

# NEK2 antibody

Catalog Number: 113102

## Product name

NEK2 antibody

## Specificity

Human, Mouse, Rat; other species not tested.

## Antibody description

NEK2 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, K-562 cells. Positive IHC detected in human testis tissue, human breast cancer tissue. Observed molecular weight by Western-blot: 48 kDa

## Preparation

This antibody was obtained by immunization of NEK2 recombinant protein (Accession Number: NM\_002497). Purification method: Antigen Affinity purified.

## Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

ELISA, WB, IHC


## Dilutions

Recommended Dilution:


WB: 1:200-1:2000

IHC: 1:20-1:200


## Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113102(NEK2 antibody) at dilution of 1:800



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:113102(NEK2 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human testis tissue slide using Catalog No:113102(NEK2 Antibody) at dilution of 1:50

# NEK2 antibody



Catalog Number: 113102

---

(under 40x lens)