

SELS antibody

Catalog Number: 115102

Product name

SELS antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

SELS Rabbit Polyclonal antibody. Positive WB detected in HeLa cells, HepG2 cells, human brain tissue. Positive IP detected in HepG2 cells. Observed molecular weight by Western-blot: 21kd

Preparation

This antibody was obtained by immunization of SELS recombinant protein (Accession Number: NM_018445). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:200-1:2000

Validations



HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:115102(SELS antibody) at dilution of 1:500



IP Result of anti-SELS (IP:Catalog No:115102, 3ug; Detection:Catalog No:115102 1:500) with HepG2 cells lysate 600ug.