

Anti-NT-proBNP (3C9) Mouse antibody



Catalog Number: 168102

Product name

Anti-NT-proBNP (3C9) Mouse antibody

Specificity

Human

Antibody description

NT-proBNP (3C9) Mouse monoclonal antibody

Preparation

Antigen: Synthetic peptide derived from human NT-proBNP.

Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300.

Storage

Store at 4°C short term. Avoid freeze / thaw cycle.

Clonality

Monoclonal


Ig Type

Mouse IgG1


Applications

ELISA

Validations



Anti-FITC mouse mAb was used to capture FITC-labeled NT-proBNP Mouse mAb (168102), HRP-labeled NGAL Mouse mAb (201289-1B5) was used for detecting.



sELISA standard Curve for NT-proBNP: Capture Antibody Mouse mAb (168102) to NT-proBNP and Detector Antibody Mouse mAb(201289-1B5) to NT-proBNP.