

# Anti-NT-proBNP (3C9) Mouse antibody



Catalog Number: 168102

## Product name

Anti-NT-proBNP (3C9) Mouse antibody

## Specificity

Human

## Antibody description

NT-proBNP (3C9) Mouse monoclonal antibody

## Preparation

Antigen: Synthetic peptide derived from human NT-proBNP.

## Formulation

Purified mouse monoclonal antibody in PBS(pH 7.4) containing with 0.03% Proclin300.

## Storage

Store at 4°C short term. Avoid freeze / thaw cycle.

## Clonality

Monoclonal


## Ig Type

Mouse IgG1


## Applications

## ELISA

### Validations



Anti-FITC mouse mAb was used to capture FITC-labeled NT-proBNP Mouse mAb (168102), HRP-labeled NGAL Mouse mAb (201289-1B5) was used for detecting.



sELISA standard Curve for NT-proBNP: Capture Antibody Mouse mAb (168102) to NT-proBNP and Detector Antibody Mouse mAb(201289-1B5) to NT-proBNP.