

ANO3 antibody

Catalog Number: 108102

Product name

ANO3 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

ANO3 Rabbit Polyclonal antibody. Positive IP detected in mouse testis tissue. Positive WB detected in mouse ovary tissue, mouse skin tissue. Observed molecular weight by Western-blot: 115kd

Preparation

This antibody was obtained by immunization of Peptide (Accession Number: NM_031418). Purification method: Antigen affinity purified.

Formulation

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IP


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IP: 1:500-1:5000

Validations



mouse ovary tissue were subjected to SDS PAGE followed by western blot with Catalog No:108102(ANO3 antibody) at dilution of 1:500



IP Result of anti-ANO3 (IP:Catalog No:108102, 4ug; Detection:Catalog No:108102 1:1000) with mouse testis tissue lysate 4000ug.