

# Human UBE2K/HIP2/LIG (GST tag) recombinant protein



Catalog Number: 519102

## General Information

### Gene Name Synonym

Ubiquitin-Conjugating Enzyme E2 K, Huntingtin-Interacting Protein 2, HIP-2, Ubiquitin Carrier Protein, Ubiquitin-Conjugating Enzyme E2-25 kDa, Ubiquitin-Conjugating Enzyme E2(25K), Ubiquitin-Conjugating Enzyme E2-25K, Ubiquitin-Protein Ligase, UBE2K, HIP2, LIG

### Protein Construction

A DNA sequence encoding the Human UBE2K/HIP2/LIG Met1-Asn200 is expressed with a GST tag at the N-terminus.

### Organism

Human

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per  $\mu\text{g}$  protein as determined by the LAL method.

### Stability

Store for up to 6 months at  $-20^{\circ}\text{C}$

### Formulation


Supplied as a  $0.2\ \mu\text{m}$  filtered solution of 50mM HEPES, 150mM NaCl, 2mM DTT, 10% Glycerol, pH 7.5.

## Usage Guide

### Storage

Store at  $<-20^{\circ}\text{C}$ , minimize freeze-thaw cycles.

### SDS-PAGE



Human UBE2K/HIP2/LIG (GST tag) recombinant protein