

U2AF1 Polyclonal Antibody

Catalog Number: 160502

Product name

U2AF1 Polyclonal Antibody

Specificity

Human, Mouse, Rat

Antibody description

Polyclonal antibody to U2AF1

Preparation

Antigen: Recombinant fusion protein containing a sequence corresponding to amino acids 1-240 of human U2AF1 (NP_006749.1).

Formulation

PBS with 0.02% sodium azide, 50% glycerol, pH7.3.

Storage

Store at -20°C. Avoid freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB, IHC, IF


Dilutions

WB 1:500 - 1:2000

IHC 1:50 - 1:100


IF 1:50 - 1:100

Validations




Western blot - U2AF1 Polyclonal Antibody

Western blot analysis of extracts of various cell lines, using U2AF1 antibody at 1:1000 dilution. Secondary antibody: HRP Goat Anti-Rabbit IgG (H+L) at 1:10000 dilution. Lysates/proteins: 25ug per lane. Blocking buffer: 3% nonfat dry milk in TBST. Detection: ECL Basic Kit. Exposure time: 90s.



Immunohistochemistry - U2AF1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded rat Intestine using U2AF1 antibody at dilution of



1:200 (40x lens).


U2AF1 Polyclonal Antibody



Catalog Number: 160502

Immunohistochemistry - U2AF1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded human oophoroma using U2AF1 antibody at dilution of 1:200 (40x lens).



Immunohistochemistry - U2AF1 Polyclonal Antibody

Immunohistochemistry of paraffin-embedded mouse liver using U2AF1 antibody at dilution of 1:200 (40x lens).