

# PALB2 antibody

Catalog Number: 113502

---

**Product name**

PALB2 antibody

**Specificity**

Human, Mouse, Rat; other species not tested.

**Antibody description**

PALB2 Rabbit Polyclonal antibody. Positive WB detected in HeLa cells. Observed molecular weight by Western-blot: 130kd

**Preparation**

This antibody was obtained by immunization of PALB2 recombinant protein (Accession Number: NM\_024675). Purification method: Antigen affinity purified.

**Formulation**

PBS with 0.02% sodium azide and 50% glycerol pH 7.3.

**Storage**

Store at -20°C. DO NOT ALIQUOT

**Clonality**

Polyclonal

**Ig Type**

Rabbit IgG


**Applications**

ELISA, WB

**Dilutions**

Recommended Dilution:

WB: 1:200-1:2000

**Validations**

HeLa cells were subjected to SDS PAGE followed by western blot with Catalog No:113502(PALB2 antibody) at dilution of 1:500