

Anti-TNFAIP8L3 antibody



Catalog Number: 176502

Product name

Anti-TNFAIP8L3 antibody

Specificity

Human, Rat, Mouse

Antibody description

Rabbit polyclonal antibody to TNFAIP8L3

Preparation

This antigen of this antibody was synthesized peptide derived from human tnfaip8l3. at aa range: 15-45

Formulation

Liquid, Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.

Storage

Store at +4°C after thawing. Aliquot store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

WB

Dilutions

WB: 1:500-1:2000

Validations




Fig1:; Western Blot analysis of 1,293T 2, mouse-brain cells using primary antibody diluted at 1:1000(4°C overnight). Secondary antibody: Goat Anti-rabbit IgG IRDye 800(diluted at 1:5000, 25°C, 1 hour)