

GCET2 antibody



Catalog Number: 110902

Product name

GCET2 antibody

Specificity

Human; other species not tested.

Antibody description

GCET2 Rabbit Polyclonal antibody. Positive WB detected in Raji cells. Positive IHC detected in human tonsillitis tissue, human lymphoma tissue. Observed molecular weight by Western-blot: 25 kDa

Preparation

This antibody was obtained by immunization of GCET2 recombinant protein (Accession Number: NM_152785). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, WB, IHC


Dilutions

Recommended Dilution:


WB: 1:200-1:1000

IHC: 1:20-1:200


Validations



Raji cells were subjected to SDS PAGE followed by western blot with Catalog No:110902(GCET2 Antibody) at dilution of 1:300



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog No:110902(GCET2 Antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemistry of paraffin-embedded human tonsillitis tissue slide using Catalog

GCET2 antibody



Catalog Number: 110902

No:110902(GCET2 Antibody) at dilution of 1:50 (under 40x lens)