

TCEB3 antibody

Catalog Number: 115902

Product name

TCEB3 antibody

Specificity

Human, Mouse, Rat; other species not tested.

Antibody description

TCEB3 Rabbit Polyclonal antibody. Positive IF detected in Hela cells. Positive WB detected in mouse liver tissue, mouse brain tissue. Observed molecular weight by Western-blot: 110kd

Preparation

This antibody was obtained by immunization of TCEB3 recombinant protein (Accession Number: NM_003198). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications

ELISA, IF, WB


Dilutions

Recommended Dilution:


WB: 1:500-1:5000

IF: 1:10-1:100

Validations



mouse liver tissue were subjected to SDS PAGE followed by western blot with Catalog No:115902(TCEB3 antibody) at dilution of 1:1000



Immunofluorescent analysis of Hela cells, using TCEB3 antibody Catalog No:115902 at 1:25 dilution and Rhodamine-labeled goat anti-rabbit IgG (red).