

# Anti-C7 Rabbit antibody



Catalog Number: 166902

## Product name

Anti-C7 Rabbit antibody

## Specificity

Human, Mouse, Rat

## Antibody description

C7 Rabbit polyclonal antibody

## Preparation

Antigen: A synthetic peptide of human C7

## Formulation

Supplied in 50nM Tris-Glycine(pH 7.4), 0.15M NaCl, 40%Glycerol, 0.01% sodium azide and 0.05% BSA.

## Storage

Store at -20°C. Stable for 12 months from date of receipt.

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications

WB, IHC-P, FC, IP

## Dilutions


WB: 1/1000-1/5000

IHC: 1/50-1/200

FC: 1/20

IP: 1/20-1/100

## Validations



Western blot detection of C7 in HeLa,293 cell lysates using C7 Rabbit pAb(1:1000 diluted).Predicted band size:94KDa.Observed band size:94KDa.