

# DFNA5 antibody



Catalog Number: 109902

## Product name

DFNA5 antibody

## Specificity

Human, Mouse; other species not tested.

## Antibody description

DFNA5 Rabbit Polyclonal antibody. Positive WB detected in mouse lung tissue, HeLa cells, Y79 cells. Observed molecular weight by Western-blot: 55 kDa

## Preparation

This antibody was obtained by immunization of DFNA5 recombinant protein (Accession Number: BC019689). Purification method: Antigen affinity purified.

## Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

## Storage

Store at -20°C. DO NOT ALIQUOT

## Clonality

Polyclonal

## Ig Type

Rabbit IgG

## Applications


ELISA, WB

## Dilutions

Recommended Dilution:

WB: 1:500-1:5000

## Validations



mouse lung tissue were subjected to SDS PAGE followed by western blot with Catalog No:109902(DFNA5 antibody) at dilution of 1:500