

# Mouse Fgf9/Fgf-9 (His tag) recombinant protein



Catalog Number: 519902

## General Information

### Gene Name Synonym

Fibroblast growth factor 9;FGF-9;Glia-activating factor;GAF;heparin-binding growth factor-9;HBGF-9;Fgf9;Fgf-9

### Protein Construction

A DNA sequence encoding the Mouse Fgf9/Fgf-9 Met1Ser208 is expressed with a 6his tag at the C-terminus.

### Organism

Mouse

### Expression Host

E. coli

## QC Testing

### Purity

>95% as determined by SDS-PAGE.

### Endotoxin

<1.0 EU per µg protein as determined by the LAL

method.

### Stability

Store for up to 6 months at -20°C

### Formulation


Supplied as a 0.2 µm filtered solution of 20mM Tris,150mM NaCl,5%Trehalose,1mM EDTA,20%glycerol□1mM DTT,pH8.5 .

## Usage Guide

### Storage

Store at <-20°C, minimize freeze-thaw cycles.

### SDS-PAGE



Mouse Fgf9/Fgf-9 (His tag) recombinant protein