

GREM2 antibody



Catalog Number: 111202

Product name

GREM2 antibody

Specificity

Human, Mouse; other species not tested.

Antibody description

GREM2 Rabbit Polyclonal antibody. Positive IHC detected in human kidney tissue, human ovary tissue.

Preparation

This antibody was obtained by immunization of GREM2 recombinant protein (Accession Number: XM_011544249). Purification method: Antigen affinity purified.

Formulation

PBS with 0.1% sodium azide and 50% glycerol pH 7.3.

Storage

Store at -20°C. DO NOT ALIQUOT

Clonality

Polyclonal

Ig Type

Rabbit IgG

Applications


ELISA, IHC

Dilutions


Recommended Dilution:

IHC: 1:20-1:200

Validations



Immunohistochemical of paraffin-embedded human kidney using Catalog No:111202(GREM2 antibody) at dilution of 1:50 (under 10x lens)



Immunohistochemical of paraffin-embedded human kidney using Catalog No:111202(GREM2 antibody) at dilution of 1:50 (under 40x lens)